

Senza

Sheath Introducer Kit



***Feel Smoothness
Throughout the Procedure***


Senza

Sheath Introducer Kit



Superior flexibility

Special material and design gives flexibility to move on guide wire and take shape in tortuous vessels



Ultra thin wall with good radial strength will allow for sheaths with smaller outer diameter but with the requisite kink resistance.




Unique valve Design

*Transparent sheath
for blood back flow visibility*


*Ultra Thin wall will reduce
post procedure hemostasis risk*

Smooth Transition

Smooth transition between dilator and sheath helps in atraumatic entry into blood vessel.



Hemostatic valve smoothness comparison



Ordering Information

Part Number*	French Size	Sheath Length (cm)	Needle Type	Guidewire Dia & material	Guidewire Tip
IMD-5F-R821	5	8	20-Gauge JELCO	0.021 inc (steel wire)	Straight

IMD-5F-R21	5	11	20-Gauge JELCO	0.021 inc (steel wire)	Straight
------------	---	----	----------------	------------------------	----------

IMD-5F-F35	5	25	18-Gauge JELCO	0.035 inc (steel wire)	J-Tip
------------	---	----	----------------	------------------------	-------

Part Number*	French Size	Sheath Length (cm)	Needle Type	Guidewire Type	Guidewire Length (cm)
IMD-6F-R821	6	8	20-Gauge JELCO	0.021 inc (steel wire)	Straight

IMD-6F-F35	6	11	18-Gauge JELCO	0.035 inc (steel wire)	J-Tip
------------	---	----	----------------	------------------------	-------

IMD-6F-F35	6	25	18-Gauge JELCO	0.035 inc (steel wire)	J-Tip
------------	---	----	----------------	------------------------	-------

Guidewire length available : 45 cm
Custom guidewires : available
Applicable for transradial and transfemoral approaches



incuMed LLP

Plot No:60, Phase-III, IDA,Cherlapally, Hyderabad, Telangana-500051.

Phone : 9000017006, 9000933920.

E-Mail : info@incumed.in

Web : www.incumed.in